Coolview®The Brightest LED Examination Lights





Coolview® Accurate Illumination For All Your Examination Needs

"Brandon Medical is a UK designer and manufacturer of world-leading equipment for operating theatres, critical care areas and primary care institutions. For over 70 years, Brandon Medical has provided healthcare professionals with reliable, high-quality and affordable medical equipment. Continuous investment in research and development of our products allows us to meet the everchanging needs of our customers, the healthcare industry and the patients it serves."



Lights to suit all of your clinical needs

Brandon Medical provides a comprehensive range of LED medical lights to suit all of your clinical needs. Brandon Medical's LED medical lights use the latest innovations in LED lighting technology to provide you with better performance and lower lifetime costs. The choice of medical lights we offer vary in light intensity, light colour, mounting structures and price so you can choose a light that meets all of your clinical examination needs.





Coolview[®] CLED50FX[™]



New Generation Examination Light

Coolview CLED50FX examination lights are the perfect light source for a wide range of procedures. Providing a high intensity, bright natural light, Coolview CLED50FX is perfect for medical consultation in primary care and medical clinics.

Powerful Light Intensity

Coolview CLED50FX emits 50,000 Lux of bright natural light at 0.5m providing an intense, clear light source for use in general consultation areas.

Perfect Colour

The colour corrected light source with a colour temperature of 4,400°K allows for accurate representation of tissue during examination.

Coolview CLED50FX has a high R_a of ≥ 93 which describes the exceptionally high quality of white light whilst the high measurement of R_g of ≥ 90 describes how accurately it illuminates tissues.

Infection Control

The fully enclosed light head is sealed to protect from dust and light liquid ingress.

Spring Balanced Arm


The FX spring balanced arm with base rotation of 360° ensures smooth, easy positioning and the arm will stay fixed into position without drifting. The modern concept is designed ergonomically to suit daily use by clinical staff.

Technical Data					
Light Intensity (at 0.5m)	50,000 Lux				
Colour Temperature	4,400°K				
Illuminated Area (at 0.5m)	180mm				
Colour Rendition Index (Ra)	≥93				
Colour Rendition Index Red Colour (Rg)	≥90				
Head Rotation	300°				
Head Angle adjustment (lamp head)	180°				
Degree of Protection	IP20				





Coolview® CLED23™ Premium Performance



areas such as neonatal and ICU departments.

Highest Light Quality

CLED23 emits 60,000 Lux of intense bright light at 0.5m, providing a pool of clear light during procedures. The light intensity can be altered through 4 easily selectable settings providing a suitable level of intensity for all procedures.

Large Illuminated Area

A large illuminated area of 210mm diameter makes it easy to examine large areas of tissue and makes it easier to target the light on the patient.

High Colour Rendition

CLED23 reproduces all of the visible colours naturally and accurately. The high R_a of ≥ 95 describes the exceptionally high quality of white light whilst the high R9 of ≥ 90 measurement describes how it accurately reproduces visible red.

Colour Control Setting

CLED23 allows the user to adjust the colour temperature of the light source. This enhances the clinician's ability to recognise and distinguish tissue conditions. The warm white (3,500°K) setting is ideal for dermatological examinations and the daylight white setting (4,500°K) is ideal for the recognition of blood vessels.








The Science Behind the Sanitizing

inhibiting the growth of micro organisms including

MRSA. This is more than an external coating and allows the lamp to be permanently resistant to contamination.



The anti-microbial substance is integrated into the light head during manufacturing. It inhibits the growth of and proliferation of microbes and bacteria, whilst the silver ions destabilise the cell membrane. This prevents the bacteria from multiplying. Tests have shown that 99% of bacteria and germs are destroyed using the Sanitized® active substance.





Environmentally Friendly

16 powerful and long lasting LEDs (30,000 hours) provide a maintenance-free light source whilst using minimal power consumption. The LEDs do not contain mercury, lead, cadmium, CFCs, POPs, VOCs or halogens.





Coolview® CLED23™ Premium Performance



Full Rotation of Light Head

The Coolview® CLED23 light head boasts 360° rotation. This allows the user to position the light beam to suit any procedure, without any restrictive stops or twisted cables. This also means the examination light can be used to provide ambient room lighting by reflecting light from the ceiling. This is a feature often used in neonatal care to provide night lighting.

Full Rotation at Ceiling Pole

Ceiling mounted models (CLED23T1C) also have unrestricted 360° rotation at the ceiling pole providing even more positioning flexibility.

Long Arm Reach & Large Angular Movements

The exceptionally long, 1200mm arm reach of the CLED23 ensures the powerful beam can be positioned all around the patient couch. This allows for 'head to toe' examination of the patient without the need for patient repositioning during examination.

Spring-Balanced Positioning

Innovative spring-balanced arm system ensures the light is positioned effortlessly, remains fixed in place and is maintenance-free.

Fully Enclosed Light Head

The light head is sealed to prevent the ingress of dust and moisture to IP20.







Intuitive Control Panel

The easy to use and ergonomic control pod allows all lighting to be adjusted with a simple press of the button.

Technical Data			
Light Intensity (at 0.5m)	60,000 Lux		
Illumination Field	210mm @ 0.5m		
LED Life	30,000 hours		
Colour Temperature	3,500°K - 4,500°K (2 steps)		
Colour Rendition Index (Ra)	≥95		
Colour Rendition Index Red Colour (Rg)	≥90		
Power Consumption	True Power 21VA Apparent Power 25VA		
IP Rating	IP20		
Head Rotation	360°		
Safety Standards	EN 60601-1, EN 60601-2-41		
Mounting Options	Ceiling, back-plate, desk, mobile, wall & rail		

^{*}Specifications may vary by +/-10%

Part Numbers	
CLED23T1C	Ceiling mounted examination light
CLED23TXBP	Back-plate mounted examination light
CLED23TXDP	Desk mounted examination light
CLED23TXMP	Mobile mounted examination light
CLED23TXWP	Wall mounted examination light
CLED23TXRP	Rail mounted examination light
CAR14443	Coolview® rail mount accessory kit

All models are supplied as standard with a 1.8m cable terminating in a UK plug.

Additional Wiring Options				
YA70439	2 metre Lead Type C Europlug to C7 Connection			
YA70438	2 metre Lead Type A Nema US 2 pin C7 Connection			



Coolview® LED Examination Lights

Mounting Options



Mounting Options

Brandon Medical's LED examination lamps are available in a variety of mountings to suit all clinical environments. If you are unsure which mounting option is most suitable, please use the table below as a guide to finding the most appropriate configuration.

Mobile

- Ideal for getting as close to the patient as possible.
- Easy to move and position.
- Use the light in multiple rooms.

Rail

- Added flexibility over a wall mounted light.
- Allows the light to be positioned nearer to the point of examination.
- Allows for easy movement along the couch length.

Ceiling

- Can be mounted anywhere in the room.
- Can be positioned out of the way.
- Ideal when space is limited.
- Longer arm reach.












Wall

- Lowest cost option.
- Can be hard wired into the wall, freeing up plugs and floor space.
- More secure than a mobile light.

Back-plate

- Should be used on non load-bearing walls for example partition walls.
- Can be hard wired into the wall, freeing up plugs and floor space.





Onesiality	Dell Manustari	Mall Manuskard	Makila Dana	Orillian Manustad	Darle Olamon
Speciality	Rail Mounted	Wall Mounted	Mobile Base	Ceiling Mounted	Desk Clamp
A&E Consulting	All Applicable	All Applicable	All Applicable	X	CLED30
A&E Minors	CLED23	CLED23	CLED23	CLED23	X
Anaesthetics	CLED30	CLED30	CLED30	CLED23	CLED30
Critical Care	CLED23/CLED10	X	CLED23	X	X
Dermatology	CLED30SV CLED23	CLED30SV CLED23	X	X	CLED30SV
	CLED10	CLED10			
Endoscopy	X	Х	CLED23	X	Х
General Practitioner	All Applicable	All Applicable	All Applicable	All Applicable	All Applicable
General Ward	CLED30 CLED50 CLED10	All Applicable	CLED30 CLED50	X	X
Gynaecology	CLED30SV CLED23 CLED10	CLED30SV CLED23 CLED10	CLED30SV	CLED23 X X	X X X
Intensive Care Unit	CLED50 CLED23 CLED10	X	CLED50 CLED23 CLED10	X	X
Neonatal Care	CLED50 CLED30SV CLED23 CLED10	X	CLED50 CLED30SV CLED23 CLED10	X	X
Obstetrics	CLED30SH CLED23	CLED30SH CLED23	CLED30SH	CLED23	X
Patient Room	CLED50 CLED30SV	CLED50 CLED30SV CLED10	X	X	X

Brandon Equipment Package

For Operating Theatres and Critical Care Areas







Medical Architectural

- Ultra clean ventilation system.
- Medical gas pendants and accessories.
- Operating tables and accessories.

Medical Audio Video

- Integrated OR solutions.
- Operating theatre 4K video.
- Digital audio video solutions.

Medical Control and Power

- Medical IT (IPS) and UPS systems.
- Intelligent theatre control panels (ITCP).

Medical Lighting

Operating theatre lights.

Exceeding Standards

Brandon Medical is certified as an ISO 13485 and ISO 9001 manufacturer, the global standard for medical device manufacturers.

Brandon Medical Coolview® lights conform to IEC 60601-1, IEC 60601-1-2 and IEC 60601-2-41





















Operating Theatre Lights and Battery Back Up
Minor Surgical Lights
Examination Lights
Theatre Control Panels
PACS Workstations
IPS Isolated Power Systems
UPS Uninterruptible Power Systems
Medical Gas Pendants and Accessories
Ultra Clean Ventilation Canopies
Digital Audio Video Solutions

NSW/ACT: 02 9723 3099 QLD: 07 3879 5999

VIC/TAS/SA/NT: 03 9532 3168 sales@salcommercial.net.au

